

STRAUMANN SIRIOS™

Capture Excellence.



STRAUMANN SIRIOS™ **CAPTURE** **EXCELLENCE**


Straumann SIRIOS™ goes beyond image capturing, redefining the essence of how you interact with dentistry and opening the door to excellence, with your patient at the center.

The synergy of our performance-driven hardware with our intuitive software is engineered to deliver peak performance across your workflow with precision and efficiency.






**REVOLUTIONARY
SCANNING PERFORMANCE**



**DESIGNED WITH YOUR
WORKFLOW IN MIND**



**UNLOCK ULTIMATE
ACCURACY**



**ELEVATED STRAUMANN®
EXPERIENCE**

“The Straumann SIRIOS™ gives me the scanning performance I need, while allowing me to quickly incorporate it into my existing workflows. Its wireless design, with optimal size, weight, and handling – makes scanning fast, easy and comfortable for my patients, whilst giving me the flexibility to adapt to individual treatment needs.”

Dr. Gustavo Harfagar

REVOLUTIONARY SCANNING PERFORMANCE

Redefine ease of use with our intraoral scanning solution. Straumann SIRIOS™ features a smart design empowering you to adapt to various clinical situations. Start scanning easily to deliver high-speed, results efficiently.

DESIGN

LIGHTWEIGHT AND COMPACT DESIGN

This wireless scanner combines cutting-edge technology with an ergonomic pen-grip design for user and patient comfort. Weighing only 245 g (including battery), it ensures effortless scanning.

TWO TIP SIZES

Select the tip size that best suits the clinical situation to support a premium care experience for your patients.

CONNECT TO THE STRAUMANN ECOSYSTEM

Connect your Straumann SIRIOS™ to the AXS platform and enjoy unlimited cloud storage and connections to Smile in a Box, your local lab and UNIQ.





EFFICIENCY

WIRELESS CHARGING

Charging is as easy as placing the device back into its holder. With Instant Wake-up, there is no waiting time, it activates the moment it is removed from the holder.

UP TO TWO HOURS OF CONTINUOUS SCANNING

Enjoy two hours of continuous scanning on a single charge, providing power and helping you stay efficient throughout your workday.

HIGH-SPEED SCANNING

Complete a full-arch scan in under 20 seconds with high-quality data capture and real-time visualization. Experience ultra-fast, smooth connectivity thanks to our advanced image compression technology.

SIMPLICITY

CALIBRATION FREE

Eliminate the need for costly calibration tools and avoid disruptions to your routine. We design the components maintain their integrity, making the device ready for precise scanning from day one and reliable over time.

DUO SMART SCANNING BUTTON

Two strategically placed scan buttons on either side of the scanner allow for seamless rotation and continuous scanning without the need to remember button positions or reassemble the tip.

EASY TO INSTALL AND MASTER

Wireless, our solution offers a quick, true plug-and-play setup, allowing you to start scanning immediately. Start scanning effortlessly with our optimized scanning technology and user-friendly software.

DESIGNED WITH YOUR WORKFLOW IN MIND

Straumann SIRIOS™ is designed to unlock your workflow to achieve excellence. Optimized for Straumann® solutions, along with integration into the new Straumann AXS® portal, it provides you with seamless and tailored workflows to respond to your needs.

Experience peace of mind with a secure environment that centralizes everything in one place, while also enabling the flexibility to work in an open-system manner.

OPTIMIZED FOR STRAUMANN® SOLUTIONS

Straumann SIRIOS™ is integrated into the Straumann® digital ecosystem with seamless workflow solutions.

IMPLANT WORKFLOWS

Streamlined scanbody matching and an intuitive interface to guide you through the workflow journey.

SURGICAL PLANNING WORKFLOW

Easily plan implant cases and design surgical guides. Export IO scans directly to coDiagnostiX® on your laptop.

RESTORATIVE WORKFLOW

Direct NOVA CAD software integration. Design temporary and final restorations and all types of models easily. Export IO scans directly to NOVA on your laptop.

ORTHODONTIC WORKFLOW

Treat aligner, revision and retainer cases from the Straumann SIRIOS™ software. Have direct patient-case connectivity to Dr Portal from ClearCorrect®.



STREAMLINED WORKFLOW THROUGH STRAUMANN® AXS

Have access to everything in one place through Straumann® AXS. Our end-to-end digital workflows are now totally seamless.

Straumann® AXS integrates software, systems and devices to streamline dental workflows. In combination with Straumann's proven implantology treatment solutions and knowhow we provide a holistic, patient-centric offering that enhances dental practice efficiency and quality of care.

SMILE IN A BOX®

Outsource guided surgery cases easily through AXS. Link IO scans to the Smile in a Box® case order form in Straumann® AXS.

SEND TO LAB WORKFLOW

Streamline your practice and easily share cases with your lab with only one click.



Discover more at:
www.straumann.com/axs

UNLOCK ULTIMATE ACCURACY

Empower your clinical decision-making with high predictability, thanks to high-resolution scanning and accurate 3D files. Experience advanced scanning technology, select your optimized workflow, and unlock numerous treatment possibilities to support delivery of an enhanced patient experience with Straumann SIRIOS™

REWORK MANAGEMENT


You can rely on scanbody recognition, validation, and streamlining operations and communications between CAD software and lab partners.

SAVE VALUABLE TIME

Straumann SIRIOS™ digital impressions can be stored indefinitely and easily reused for future restorations. This may eliminate errors and the need for new scans saving valuable time and reducing rework costs.

CLINICAL OUTCOMES

Enables accurate scanning empowering you to achieve optimal clinical outcomes. Enhanced fit and function, compared to traditional impressions, combined with the ability to offer a broad range of treatment options, intended to support patient satisfaction.



Engage with your patients, track their records and show treatment options directly on screen.



Smart Scanning Button

Start scanning with a single touch and monitor your progress in real-time with dynamic lighting guidance.



Mode Switch Button

Easily navigate through different workflow stages with a click.



Status Indicator

Check scanning status and progress with lighting indication.



Battery Indicator

Be reminded about battery level status with lighting indication.



THE ELEVATED STRAUMANN® EXPERIENCE

From Straumann®, you can expect nothing less than premium service and customer experience that follow you since the beginning of your journey with us.

SUPERIOR TRAINING AND EDUCATION

Our commitment to education is a cornerstone of effective patient care and practice success. We offer a wide range of educational programs, including classroom courses, eLearning, webinars – all integrated into one digital platform – worldwide.

For Straumann SIRIOS™ you can also count on a unique e-Learning Masterclass onboarding. It features dynamic and innovative content to provide you with an engaging learning experience.

SUPPORT YOU CAN RELY ON

Have access to a customer support program that offers quick issue resolution, availability in multiple geographic locations, and local language assistance.

Benefit from a comprehensive customer care package that includes extended warranty, full accidental coverage, and all the valuable Straumann® customer care services you already know and trust.



Straumann® is the flagship brand of the world's largest implant manufacturer group – the Straumann Group. Renowned as one of the most iconic brands in dentistry, we have a long-standing legacy of pioneering significant advancements in the field. We are at the forefront of expanding the scope of dental care by relentlessly and consistently providing excellence and transformational innovation to a global community of forward-thinking dental professionals.



STRAUMANN SIRIOS™ CAPTURE **EXCELLENCE.**



NAMLIT.1752 1/25 V1 PMR

Straumann North American Headquarters

Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Phone 800/448 8168 (US) · 800/363 4024 (CA)
Fax 978/747 2490
www.straumann.us · www.straumann.ca

ifu.straumann.com
<https://ifu.allied-star.com/IOS>

© Straumann USA, LLC 2025. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.

Note: For recognized legal manufacturer, refer to the product label.

